

Preventative Maintenance Guide Dansensor CheckPoint 3 (S)

Technical Support: 01223 492087 techsupport@tendringpacific.com

Dansensor Preventative Maintenance Guides help you to care more effectively for your equipment during the time between planned services.

For Safety Instructions, Trouble Shooting & Operation - please refer to the User Guide.

Cleaning and Maintenance

In addition to the Annual Service & Calibration, it is recommended to carry out the following maintenance checks at the interval times stated below.

During use of the device, dust, debris, and moisture may be sucked through the measuring system.

Every Day:


Water Trap Filter and Needle - Check and Replace (as required) The water trap filter can become blocked due to debris, or torn by the debris, because of internal moisture. It will need to be replaced regularly.

The needle may sometimes touch the product being measured, or the needle may suck fluid from the packaging from which the measurement is taken. If this happens, it will be necessary to replace both the needle and the filter.

Once a Week:

Sample Hose - Check and Replace (as required)

If the needle and water trap filter have been replaced and a blockage error persists, then replace the sample hose system. To replace the sample hose, push downwards on the sleeve (4) to release the hose (3), then pull the hose upwards and out of the connector. Insert the new hose by pushing it into the connector. Check for proper connection, by pulling the hose slightly upwards. 20.9% O₂ air calibrations are not performed on the Dansensor CheckPoint 3 (S) model.



If your CheckPoint 3 (S) is showing readings out of specification, please call technical support for assistance.

The unit may need repairing or recalibration.

Dansensor CheckPoint 3 Consumable Parts

Consumable Kit, comprises of: Part No: 360408

1 x Septum, 15mm grey soft (1000pcs)

1 x Needle kit (10pcs)

1 x Sampling kit

1 x Filter water trap (10pcs)

Filter water trap kit (10pcs) Part No: 310335

Filter water trap kit (100pcs)

Needle kit (10pcs)

Part No: 310339

Part No: 340532

Needle kit stainless steel (12pcs)

Part No: 340566

Sampling kit

Part No: 340516

Septums, grey (1000pcs)

Part No: 310337

Septums, grey (10000pcs)

Part No: 220157

To order consumables please contact Ann Walker on 01223 495484 or email ann.walker@tendringpacific.com

